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
Mathews, C.K m.fl. Biochemistry. Benjamin/Cummings. 1997

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

Fig. 5.4. Structure of plant chlorophylls. (The arrangement of double bonds in the macrocyclic ring of chlorophyll is not fixed. Resonance of the conjugated double-bond system makes it possible to write down a number of mesomeric structures, one of which is shown above. The positions of the two co-ordinate bonds (→) and the two covalent bonds between the Mg²⁺ and pyrrole nitrogens, whilst being appropriate to the mesomer shown, are also not fixed for the same reason. The two positive charges on the Mg²⁺ are balanced by two negative charges which are shared amongst the four pyrrole nitrogens; thus chlorophyll is electrically neutral.)


Textual No

XVII

Phycoerythrobilin

XVIII

Phycocyanobilin

XIX


Fig. 5.6. Structures of phycobilins.


Fig. 5.5. Structures of some of the main plant carotenoids (see Table 5.1). Carotenoids VI–XI have the same central segment; Cs 7 and 7' are bonded to Cs 6 and 6' of the left-hand and right-hand rings respectively, at the same level in the figure.


Fig. 5.7. Absorption spectra of selected members of the three groups of photosynthetic pigments, the chlorophylls, the carotenoids and the phycobilins.


Fig. 5.8. Curve analyses of the red peak in the absorption spectra, taken at -196°C , of fragments of stroma lamellae (spectrum A) which constitute a chloroplast fraction enriched in PS I and of fragments of grana lamellae (spectrum B) which constitute a chloroplast fraction enriched in PS II (adapted from French, Brown and Lawrence, 1972).

Table 5.1. Distribution of the major photosynthetic pigments in the plant kingdom

Organism	Chlorophylls					Phycobiliproteins		Carotenoids		
	<i>a</i>	<i>b</i>	<i>c</i> ₁	<i>c</i> ₂	<i>d</i>	Phycoerythrin	Phycocyanin	Carotenes	Xanthophylls	
Higher plants, pteridophytes and bryophytes	+	+	-	-	-	-	-	β -carotene α -carotene	lutein violaxanthin neoxanthin	
Algae	Chlorophyceae (green)	+	+	-	-	-	-	β -carotene	lutein violaxanthin neoxanthin	
	Euglenophyceae	+	+	-	-	-	-	β -carotene	neoxanthin diadinoxanthin	
	Phaeophyceae (brown)	+	-	+	+	-	-	β -carotene	fucoxanthin violaxanthin	
	Chrysophyceae (yellow-brown)	+	-	+	+	-	-	β -carotene	fucoxanthin	
	Xanthophyceae (yellow-green)	+	-	-	-	-	-	β -carotene	neoxanthin diadinoxanthin	
	Bacillariophyceae (diatoms)	+	-	+	+	-	-	β -carotene	neoxanthin diadinoxanthin fucoxanthin	
	Cryptophyceae	+	-	-	+	-	+	+	α -carotene β -carotene	alloxanthin
	Rhodophyceae (red)	+	-	-	-	+	+++	+	α -carotene β -carotene	lutein zeaxanthin
	Cyanophyceae ^(a) (blue-green)	+	-	-	-	-	+	+++ ^(c)	β -carotene	echinenone myxoxanthophyll zeaxanthin
Prochlorophyta ^(b) (genus: Prochloron)	+	+	-	-	-	-	-	β -carotene	zeaxanthin	

+ and - = presence and absence respectively.

^(a) The taxonomic position of the blue-green algae is at present in a state of flux: some authors believe that they should be classed as bacteria because they are prokaryotic and therefore refer to them as cyanobacteria; other authors believe that they should remain within the algal classification because their mechanism of photosynthesis is more plant-like than bacteria-like (they produce oxygen and have two photosystems). Being prokaryotes they do not have chloroplasts.

^(b) Prochloron constitutes the only known genus within the new algal division Prochlorophyta. It is symbiotically associated with colonies of didemnid ascidians and lives in warm marine habitats. It is prokaryotic but its photosynthetic mechanism is plant-like in that it evolves oxygen and has two photosystems which are built into the membranes of paired thylakoids.

^(c) Allophycocyanin is also present in small quantities.